

Revised: 07–June–2007

Leadmium[™] Green AM Dye for Intracellular Detection of Lead and Cadmium





Introduction

Heavy metals are ubiquitous environmental contaminants and exhibit widespread toxicity. Although standard assays exist for metal ions in solution, specific indicators to detect intracellular levels of lead or cadmium were not previously available. The Intracellular Lead and Cadmium Detection Kit employs the Leadmium[™] Green AM dye as a specific indicator of lead or cadmium in cells. After the cells are loaded with dye, nonspecific esterases cleave the acetoxymethyl (AM) group, resulting in a charged form of the Leadmium[™] Green dye that leaks out of cells far more slowly than the parent compound. The calciuminsensitive Leadmium[™] Green dye becomes fluorescent in the presence of nanomolar levels of lead and micromolar levels of cadmium. This dye can be used in both microscopy and flow cytometry applications.

Product Information:

Components

• Leadmium[™] Green AM dye (Component A), 50 µg in each of 5 vials

Storage and Handling

Upon receipt, store dye at $\leq -20^{\circ}$ C. To minimize condensation, allow vials to warm to room temperature before opening.

Once prepared, solutions of the LeadmiumTM Green AM dye in DMSO should be used within 1 day.

Figure 1. Background fluorescence of Leadmium™ Green dye in the Jurkat cell line. Cells in saline were loaded with Leadmium™ Green AM dye and washed. The singlecolor fluorescence was collected using 488 nm excitation and a 525/10 nm bandpass filter. Because no intracellular lead or cadmium was present in the cells, this histogram indicates the background fluorescence contributed by Leadmium™ Green dye.

Spectral Characteristics

The approximate excitation and emission maxima of Leadmium[™] Green dye are 490 nm and 520 nm, respectively. Cells labeled with the Leadmium[™] Green dye can be analyzed by flow cytometry using 488 nm excitation and measuring fluorescence emission at 520 nm.




Figure 2. Fluorescence of Jurkat cells containing PbCl₂ in the Jurkat cell line. Cells in saline were loaded with LeadmiumTM Green AM dye and washed and then incubated in the presence of 1 μ M PbCl₂ in saline and 1 μ M ionomycin. The single-color fluorescence was collected using 488 nm excitation and a 525/10 nm bandpass filter. The histogram shows cells positive for intracellular lead.

Before You Begin

Materials Required but Not Provided

- Saline (0.85% NaCl);
- **Note:** phosphate-buffered saline (PBS) should not be used DMSO, high-quality, anhydrous

Optional Reagents

- Solutions of PbCl
- Solutions of CdCl
- Tetrakis-(2-pyridylmethyl) ethylenediamine (TPEN, Cat. no. T1210)
- Ethanol
- Ionomycin, calcium salt (Cat. no. I24222)
- Dead-cell stain such as propidium iodide (Cat. no. P3566, 1.0 mg/mL solution in water); SYTOX[®] Red dead cell stain for 633 or 635 nm excitation (Cat. no. S34859, 5 μM solution in DMSO); or SYTOX[®] Blue dead cell stain for flow cytometry (Cat. no. S34857, 1 mM solution in DMSO)

Reagent Preparation


1.1 Remove one vial of Leadmium[™] Green AM dye from freezer and let come to room temperature; protect from light.

1.2 Make the LeadmiumTM Green AM stock solution by adding 50 μ L DMSO to one vial of dye. Mix well and protect from light.

1.3 Prepare a Leadmium[™] Green AM dye working solution by diluting Leadmium[™] Green AM stock solution (prepared in step 1.2) in saline 1:10 (i.e., adding 20 µL stock solution dye to 180 µL saline). Protect dye solutions from light at all times.

1.4 Warm saline to 37°C. This will be used for washing the samples between various steps.

1.5 OPTIONAL: If the experiment will include comparative controls (treatment of cells with known concentrations of $PbCl_2$), a recommended starting concentration range for $PbCl_2$ is 625 nM



Leadmium Green dye fluorescence

Figure 3. Dual-color scatter plot showing two populations. Jurkat cells were loaded with LeadmiumTM Green AM dye and washed. The sample was then incubated in the presence of 1 μ M PbCl₂ (in saline) and 1 μ M ionomycin. After washing the sample, it was incubated in the presence of PI. Dual-color fluorescence was collected using 488 nm excitation and 525/10 nm and 610/10 nm bandpass filters. This results in visualization of two populations: dead cells (positive for PI) and live cells that contain lead (positive for LeadmiumTM Green dye and negative for PI).

down to 1 nM. Make 625 nM PbCl₂ by adding 250 μ L of 10 μ M PbCl₂ to 3.75 mL saline. Make 9 further serial dilutions of this to arrive at 1 nM PbCl₂. Store these solutions at 37°C until required. You will need 1 mL of each concentration for one treatment. Adjust the volumes transferred in the serial dilutions if replicates are desired.

1.6 OPTIONAL: If the experiment will include comparative controls (treatment of cells with known concentrations of CdCl₂), a recommended starting concentration range for CdCl₂ is 250 μ M down to 1 μ M. Make 250 μ M CdCl₂ by adding 1 mL of 1 mM CdCl₂ to 3 mL saline. Make 8 further serial dilutions of this to arrive at 1 μ M CdCl₂. Store these solutions at 37°C until required. You will need 1 mL of each concentration for one treatment. Adjust the volumes transferred in the serial dilutions if replicates are desired.

1.7 OPTIONAL: Make a 10 mM working solution of TPEN by adding 1.2 mL of ethanol to 5 mg TPEN. Mix well and store at room temperature.

Protocol for Flow Cytometry Testing

This protocol was optimized using Jurkat cells, a human T-cell line. To achieve optimal results, experimental parameters such as suspension buffers, media, and cell concentration should be adjusted depending on cell type and culture conditions used. This kit contains enough dye for 200 flow cytometry assays, based on a 1 ml assay volume.

Note: Phosphate buffers should be avoided due to metal salt insolubilities.

Staining Cells

2.1 Harvest the cell sample(s) to be assayed and pellet them by centrifugation.

Note: The basic experiment includes both the experimental cell sample (i.e., those cells for which the lead or cadmium status is unknown) and the comparative control cell sample (i.e., cells that will be subjected to lead or cadmium at known concentrations).

2.2 Remove supernatant and resuspend each sample pellet in warm saline to wash cells.


2.3 Collect the cells by centrifugation and resuspend the cell pellet (s) in warm saline at 1×10^6 cells per mL.

2.4 Set up the following tubes, each containing 1 mL cell suspension:

- **2.4a** A tube (or series of tubes) containing the experimental cells.
- **2.4b** A series of tubes containing the comparative control cells.

2.5 Add 4 μ L of LeadmiumTM Green AM dye working solution (prepared in step 1.3) to all tubes and vortex briefly to mix.

2.6 Incubate tubes in a water bath at 37°C for 30 minutes, protected from light. This incubation time may be extended for up to 60 minutes, if desired.



PbCl₂ concentration (nM)

Figure 4. Approximate level of intracellular lead in experimental cells may be determined by comparing MFI values obtained to MFI values from the comparative control samples. Jurkat cells were loaded with Leadmium[™] Green AM dye and washed. Samples were then incubated in the presence of PbCl₂ solutions in the range of 5–625 nM, and 1 µM ionomycin. After washing the samples, propidium iodide (PI) was added. Dual-color fluorescence was collected using 488 nm excitation and 525/10 nm and 610/10 nm bandpass filters. MFI measurements were recorded from a region on the dual-color plot drawn around the dye-positive and PI-negative population, and then plotted.

2.7 Add 4 mL of warm saline to each tube.

2.8 Pellet cells by centrifugation.

- 2.9 Set up the tubes for second incubation as follows:2.9a For experimental samples, add 1 mL warm saline to each tube.
 - 2.9b For comparative control samples, add 1 mL of each concentration of PbCl₂ or CdCl₂ (see steps 1.5 and 1.6) and 2 μL 500 μM ionomycin to this series of tubes.
 - **2.9c** For single-color dead-cell compensation tube, add 1 mL warm saline.
 - **2.9d** For single-color Leadmium[™] Green dye compensation tube add 1 mL either 1 mM PbCl₂ or 1 μM CdCl₂ and 2 μL ionomycin.
 - 2.9e For dye background tube, add 1 mL saline.

2.10 Vortex all tubes briefly and incubate in a 37°C water bath for 30 minutes, protected from light.

2.11 Add 4 mL of warm saline to each tube.

2.12 Pellet cells by centrifugation.

2.13 Add 0.5 mL saline to all tubes and mix well.

2.14 Add appropriate amount of dead-cell stain to all tubes (except single-color Leadmium[™] Green dye compensation tube) and vortex briefly to mix.

2.15 Incubate tubes for 5 minutes at room temperature, protected from light.

2.16 Put tubes on ice and protect from light. Analyze on the cytometer within one hour.

Flow Cytometric Analysis

3.1 Set the instrument for 488 nm excitation and use appropriate bandpass filters. A dual-parameter fluorescence plot with a region around dye-positive and dead-cell–negative events may be used to gate out dead cells.

3.2 Run unstained cells to get main population of cells on the FSC/SSC plot and adjust the voltages of the fluorescence plots to put unstained cells in the first decade.

3.3 Run single-color compensation tubes and adjust settings.

3.4 Run dye background tube. Single-color Leadmium[™] Green dye fluorescence will give data similar to that in Figure 1. Exclude dead cells by making a region around Leadmium[™] Green dye–positive and dead-cell stain negative cell population.

3.5 Analyze remaining tubes. Collect and record MFI of the dye-positive and dead-cell-negative region in the dual plot. Single-color fluorescence histogram of the intracellular lead is shown in Figure 2. Dual-color fluorescence of intracellular lead in live cells is shown in Figure 3. Graphs of the fluorescence values obtained when LeadmiumTM Green-loaded Jurkat cells were treated with various concentrations of either PbCl₂ or CdCl₂ are shown in Figures 4 and 5 respectively.

3.6 To confirm that the observed fluorescence is due to the presence of lead or cadmium, it is possible to chelate the heavy metals. To do this, allow the tubes come to room temperature, add 2.5 μ L TPEN, and vortex briefly to mix.

3.7 Incubate 5 minutes at room temperature.

3.8 Analyze the tubes at the same settings used for testing. **Note:** TPEN is a cell-permeant, heavy metal chelator. Samples treated with TPEN should show a drop in fluorescence.




Figure 5. Approximate level of intracellular cadmium in experimental cells may be determined by comparing MFI values obtained to MFI values from the comparative control samples. Jurkat cells were loaded with LeadmiumTM Green AM dye and washed. Samples were then incubated in the presence of $CdCl_2$ solutions in the range of 5–300 μ M, and 1 μ M ionomycin. After washing the samples, propidium iodide (PI) was added. Dual-color fluorescence was collected using 488 nm excitation and 525/10 nm and 610/10 nm bandpass filters. MFI measurements were recorded from a region on the dual-color plot drawn around the LeadmiumTM Green dye–positive and PI-negative population, and then plotted.

Product List Current prices may be obtained from our Web site or from our Customer Service Department.

Cat #	Product Name	Unit Size
A10024	Leadmium™ Green AM dye *for intracellular detection of lead and cadmium	5 × 50 µg
M36353	Measure-iT™ Lead and Cadmium Assay Kit *1000 assays*	1 kit
P3566	propidium iodide *1.0 mg/mL solution in water*	10 mL
P21493	propidium iodide *FluoroPure™ grade*	100 mg
S34859	SYTOX® Red dead cell stain *for 633 or 635 nm excitation* *5 µM solution in DMSO*	1 mL
S34857	SYTOX® Blue dead cell stain *for flow cytometry* *1000 assays* *1 mM solution in DMSO*	1 mL

Contact Information

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Please visit our website - probes.invitrogen.com - for the most up-to-date information.

Molecular Probes, Inc. 29851 Willow Creek Road, Eugene, OR 97402 Phone: (541) 465-8300 • Fax: (541) 335-0504	Invitrogen, Ltd. 3 Fountain Drive Inchinnan Business Park Paisley PA4 9RF, UK
Customer Service: 6:00 am to 4:30 pm (Pacific Time) Phone: (541) 335-0338 • Fax: (541) 335-0305 • probesorder@invitrogen.com	Phone: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260 Email: euroinfo@invitrogen.com Technical Services: eurotech@invitrogen.com
Toll-Free Ordering for USA: Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228	
Technical Service: 8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 • Toll-Free (800) 438-2209 Fax: (541) 335-0238 • probestech@invitrogen.com	Molecular Probes products are high-quality reagents and materials in-

Invitrogen European Headquarters

tended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The manufacture, use, sale or import of this product may be subject to one or more patents or pending applications owned or licensed by Invitrogen Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity) in a manner consistent with the accompanying product literature. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents that are owned or controlled by Invitrogen Corporation and/or Molecular Probes. Inc. which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2007, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.