

Poly(I:C)

High Molecular Weight

Synthetic analog of dsRNA - TLR3 ligand

Catalog # tlr1-pic, tlr1-pic-5

For research use only

Version # 11C21-MM

PRODUCT INFORMATION

Content:

Poly(I:C) is provided lyophilized and is available in two sizes:

- 10 mg: (catalog # tlr1-pic)
- 50 mg: (catalog # tlr1-pic-5)
- 10 ml or 2 x 25 ml sterile endotoxin-free physiological water (NaCl 0.9%)

Storage:

- Poly(I:C) is shipped at room temperature and should be stored at 4°C.
- Upon resuspension, prepare aliquots of Poly(I:C) and store at -20°C for long term storage. Store at 4°C for short term storage.
- Lyophilized product is stable 1 year at 4°C when properly stored. Resuspended product is stable 1 month at 4°C and 1 year at -20°C. Avoid repeated freeze-thaw cycles.

Quality control:

- Absorbance spectrum
- Gel retardation (Size: 1.5-8 kb)
- Human TLR3 (hTLR3) activity tested using HEK-Blue™ hTLR3 cells
- Endotoxin level: <0.001 EU/μg (measurement by kinetic chromogenic LAL assay)

DESCRIPTION

Polyinosinic-polycytidylic acid (poly(I:C)) is a synthetic analog of double stranded RNA (dsRNA), a molecular pattern associated with viral infection. Both natural and synthetic dsRNAs are known to induce type I interferons (IFN) and other cytokines production. Poly(I:C) is recognized by Toll-like receptor 3 (TLR3)^{1,2}. Upon poly(I:C) recognition, TLR3 activates the transcription factor interferon regulatory factor 3 (IRF3), through the adapter protein Toll-IL-1 receptor (TIR) domain-containing adapter inducing IFN-β (TRIF, also known as TICAM-1)³. Activation of IRF3 leads to the production of type I IFNs, especially IFN-β. A second pathway involves the recruitment of TNF receptor-associated factor 6 (TRAF6) or receptor interacting protein 1 (RIP1), with the subsequent activation of the transcription factors NF-κB and AP-1⁴. Activation of this pathway triggers the production of inflammatory cytokines and chemokines such as TNF-α, IL-6 and CXCL10. Poly(I:C) is also recognized by the cytosolic RNA helicases retinoic acid-inducible protein 1 (RIG-I) and melanoma differentiation-associated gene 5 (MDA-5)⁵.

1. Alexopoulou L. et al., 2001. Recognition of double-stranded RNA and activation of NFκB by Toll-like receptor 3. *Nature*, 413(6857):732-8. 2. Matsumoto M. et al., 2002. Establishment of a monoclonal antibody against human Toll-like receptor 3 that blocks double-stranded RNA-mediated signaling. *BBRC*, 293(5):1364-9. 3. Yamamoto M. et al., 2003. Role of Adapter TRIF in the MyD88-Independent Toll-Like Receptor Signaling Pathway. *Science* 301: 640. 4. Kawai T & Akira S., 2008. Toll-like receptor and RIG-I-like receptor signaling. *Ann NY Acad Sci*. 1143:1-20. Review. 5. Kato H. et al., 2006. Differential roles of MDA5 and RIG-I helicases in the recognition of RNA viruses. *Nature*. 441(7089):101-5. 6. Schindler U. & Baichwal VR., 1994. Three NF-κB binding sites in the human E-selectin gene required for maximal tumor necrosis factor alpha-induced expression. *Mol Cell Biol*, 14(9):5820-5831.

METHODS

Preparation of sterile stock solution (1 mg/ml)

- Stimulation of TLR3 can be achieved with 30 ng - 10 μg/ml Poly(I:C).
- Add 10 ml of the endotoxin-free physiological water provided to the 10 mg Poly(I:C) vial or 50 ml to the 50 mg Poly(I:C) vial to obtain a solution at 1 mg/ml.
 - Mix the solution by pipetting up and down.
 - Heat the mixture for 10 minutes at 65 - 70°C. Allow the solution to cool for 1 hour at room temperature to ensure proper annealing.

Example of *in vitro* activation of TLR3 with Poly(I:C) using HEK-Blue hTLR3 cells

Poly(I:C) can be used to stimulate hTLR3 in HEK-Blue™ hTLR3 cells. HEK-Blue™-hTLR3 cells are designed for studying the stimulation of hTLR3 by monitoring the activation of NF-κB. Stimulation with a TLR3 ligand activates NF-κB and AP-1 which induces the production of SEAP. Levels of SEAP can be easily determined with QUANTI-Blue™ (a detection medium that turns purple/blue in the presence of alkaline phosphatase). A typical stimulation curve is given figure 1 overleaf.

- Prepare a HEK-Blue™ hTLR3 cell suspension (250,000 cells/ml) in DMEM, 4.5 g/l glucose, 10% (v/v) heat-inactivated fetal bovine serum (30 min at 56°C), 50 U/ml penicillin, 50 μg/ml streptomycin, 100 μg/ml Normocin™, 2 mM L-glutamine.
- In a 96-well plate, add 180 μl of the HEK-Blue™ hTLR3 cell suspension per well.
- Stimulate cells with 30 ng -10 μg/ml Poly(I:C) for 6 to 24 hours.
- Determine poly(I:C) stimulation on TLR3 by assessing reporter gene expression using QUANTI-Blue™ or HEK-Blue™ detection.

Note: *InvivoGen provides also a low molecular weight poly(I:C), named poly(I:C)-LMW (see "Related Products"), with an average size of 0.2-1 kb that may activate the immune system differently.*

RELATED PRODUCTS

Product	Catalog Code
HEK-Blue™ hTLR3	hkb-htlr3
QUANTI-Blue™	rep-qbl
HEK-Blue™ Detection	hb-det1
Poly(I:C)-LMW	tlr1-picw
Poly(A:U)	tlr1-pau
Poly(I:C) Rhodamine	tlr1-picr

TECHNICAL SUPPORT

Toll free (US): 888-457-5873
Outside US: (+1) 858-457-5873
Europe: +33 562-71-69-39
E-mail: info@invivogen.com
Website: www.invivogen.com



3950 Sorrento Valley Blvd. Suite 100
San Diego, CA 92121 - USA




Figure 1. HEK-Blue™ hTLR3 cells were stimulated with increasing concentrations of Poly(I:C). After 18h incubation, NF-κB-induced SEAP activity was assessed using QUANTI-Blue™.

Performance of this assay was validated under optimized conditions in a 96-well plate using QUANTI-Blue™.

Poly(I:C) EC₅₀ = 70 +/- 10 ng/ml

Response Ratio = 15

Optimum cell number = 50,000 cells/well

TECHNICAL SUPPORT

Toll free (US): 888-457-5873

Outside US: (+1) 858-457-5873

Europe: +33 562-71-69-39

E-mail: info@invivogen.com

Website: www.invivogen.com



3950 Sorrento Valley Blvd. Suite 100
San Diego, CA 92121 - USA